

# **Technologically Advanced**

Combining together a stunning ultra-slim design, large touch screen, and a rich set of features, the iM3 presents itself a next generation vital signs monitor.



Portable Design



360° Alarm Light



Barcode Scanner Support



8" Touch Screen Display



No-Fan Design



USB, OTG, LAN, and DCHP



Built-In WiFi

#### Features:

120 hour trend review
1200 NIBP measurement storage
200 alarm storage
Data upload via USB to EMR system
3 quick-action keys
Customizable shortcuts
Multiple working modes

#### **Standard Parameters:**

SpO2, NIBP, PR

#### **Optional Parameters:**

SunTech NIBP, Quick temp, Infrared Ear Temp

#### **Built-in Memory**

Monitoring Mode (for each patient file)

200
Alarm
1200
NIBP Measurements
120 h
Trend Review

Spot & Round Mode

16,000,000 Measurements

# **Innovative Design**

#### **Adaptive Work Modes**

#### **Round Mode**

Designed to make daily use more efficient
Easily upload data through USB or downloaded from the EMR
Features color-coded list
Windows interface design

#### **Monitoring Mode**

Designed for continuous patient monitoring Real-time data, alarms, and history trend

#### **Spot Mode**

Designed for spot-check applications Real-time vitals

# TEMP Module Options







Infrared Ear TEMP

Quick TEMP

Covidien Quick Temp

## **Standard Accessories**

- NIBP Cuff, Adult, 27cm-35cm, reusable with NIBP tube
- Adult Finger Clip Sensor 1 meter with 2 meter extension cable
- Internal Li-Ion Battery
- Operator's Manual

## **Sensor Options**

- SPO2 Finger Sensor, Adult, 1m, reusable
- SpO2 Wrap Sensor, Neonate, 1m, reusable
- SpO2 Silicone Soft-tip Sensor, Pediatric, 1m, reusable
- SpO2 Ear Clip Sensor, Adult/Pediatric, 1m, reusable
- SpO2 Silicone Soft-tip Sensor, Adult, 1m, reusable

## For More Information

**C**: 877-842-9753

E: info@cardiotechmed.com

5963 Olivas Park Dr Suite F, Ventura, CA 93003

#### **Hours of Operation:**

MON - FRI SAT - SUN 6AM-6PM CLOSED