

Expert Medical Care

The CardioTech GT-300 EKG Machine boasts exceptional cost performance with both accurate clinical performance as well as user-friendly design.



Built-in Thermal Recorder



USB Support



Barcode Scanner Support



Shortcut Keys



8-inch LCD Screen



Alphanumeric Keyboard



PDF, SCP, FDA-XML, DICOM

Features:

On-board PDF Creator
Alphanumeric Keyboard
One-Touch Operation
Built-In WiFi
LCD Display
Rechargeable Li-on Battery
Printer
Optional PC Software


















User Friendly Design

SmartECG Viewer



Data Management

100 ECG internal memory (can be enlarged by USB flash disk)
Optional SmartECG viewer
DATA transmitting to PC via LAN/WIFI
PDF, SCP. FDA-XML and DICOM data export options
Bi-directional data communication with DMS

Exceptional Clinical Performance

Weak signal detection with 24 bits A/D converter technology
Anti-noise technology with comprehensive filters
High CMRR over 115dB
120s Real-time ECG waveforms freezing and review
Automatic measurements and interpretation tested by CSE/AHA/MIT database
Extended printing when arrhythmia is detected
Integrated advanced ECG analysis: Vector ECG

What's the difference?

The CardioTech GT-300 EKG and GT-350 EKG are suitable for clinics, hospitals, and offices. While the machines are very similar, the GT-350 does include a few enhanced features.

GT-300

Short-cut keys only 100 EKG memory

GT-350

Touch screen and short-cut keys 200 EKG memory

For More Information

C: 877-842-9753

E: info@cardiotechmed.com

5963 Olivas Park Dr Suite F, Ventura, CA 93003

Hours of Operation:

MON - FRI 6AM-6PM SAT - SUN CLOSED